


PI-80X Tissue DNA protocol

Genomic DNA isolation from Protease K digested mouse tissue




KURABO PI-80X is an automated DNA isolation system with throughput of 48 samples at a time. This data sheet shows an example of DNA isolation from protease digested mouse tissues. The specially designed 8-hole tube makes isolation process easy.



Experiment

Sample	Protease K digested mouse tail, brain, liver, kidney and heart.
Sample amount	Tail: 10mg The other tissues: 20mg
Isolation system and protocol	PI-80X Tissue DNA protocol
Reagents kit and consumables	NR-201 Sample tube, PT-8000; DNA tube, NT-8000
Chemical principle	Lysis of tissue: Protease K digestion (pre-treatment) Purification: Phenol extraction DNA isolation: Precipitation by alcohol
Process	<p>Pre-treatment: Digest with 0.2mg/ml of protease K at 50-55C for 15hr. Reaction component is as follow;</p> <p style="margin-left: 40px;">0.25ml 0.4mg/ml protease K in Protease dissolving reagent (No.2, NR2025) 0.25ml Suspending reagent (No. 10, NR10025) Total volume 0.5ml</p> <p>DNA isolation with PI-80X: Step1. Denature protein and remove debris Step2. Precipitate DNA Step3. Wash DNA Step4. Dry up Step5. Dissolve DNA in DNA dissolving reagent (KURABO, NR7025) Final volume 100ul</p>

Workflow of Tissue DNA protocol



Processing time 48 samples: 2.7hr
 (PI-80X) 48 samples: 2.1hr with optional heater unit.

Analysis methods

Yield and purity calculation Absorbance of 260nm was measured for each DNA solution by spectrometer. DNA yield was calculated as follows:
 $A_{260} \times 50 \times \text{dilution factor} \times \text{final volume}$.
 DNA purity was evaluated by the ratio of A_{260}/A_{280} .

Electrophoresis Each 5µl of DNA solution was loaded on a 0.7% agarose gel.

Restriction enzyme digestion Each 1ug of DNA isolated from mouse tail was digested with 10units of enzymes for 15hours at the optimal temperature, respectively.
 Restriction enzyme: *EcoRI*, *KpnI*, *SacI* and *BamHI*


PCR amplification Template: 200ng of isolated DNA from mouse tail
 Target gene: Mouse ICAM 1, exon 6-7, amplicon size 535bp
 Polymerase: AmpliTaq Gold® DNA Polymerase (0.6U)
 PCR condition: 94C, 5min x 1 cycle
 94C, 20sec/68C, 40sec/72C 1.5min x 35 cycles
 72C, 7min x 1 cycle
 Reaction volume: 25ul
 Electrophoresis: 10ul of 25ul reaction mixture was loaded on a 1.8% agarose gel.

Results


Yield and purity

Sample	Tail	Brain	Liver	Kidney	Heart
Yield (ug)	20.1/10mg	30.8/20mg	22.7/20mg	21.5/20mg	17.5/20mg
Purity (A260/A280)	1.87	1.93	1.70	1.93	1.89

Electrophoresis



Each 5ul of 100ul DNAsolution isolated from mouse tissues.
 The starting amount is;
 tail: 10mg
 the other tissues: 20mg

Restriction enzyme digestion


Each 1 μ g of DNA isolated from mouse tail was digested with restriction enzymes.

M : λ HindIII size maker


C : undigest DNA

1 : *EcoRI*

2 : *KpnI*

3 : *Sac I*

4 : *BamHI*

PCR amplification


M: ϕ X174 HinII size marker

C+: Positive control

C-: Negative control

1-7: Amplicon of isolated DNA

Ordering information

DNA isolation system PI-80X
Tissue DNA protocol

Reagents kit Tissue DNA isolation reagents kit NR-201 (approx. for 650 preps)

Reagent name	Reagent No. (Code)	Content
Proteinase dissolving reagent	2 (NR2025)	1
Deproteinization reagent	3 (NR3025)	1
Precipitating reagent	4 (NR4050)	1
Washing reagent	5 (NR5050)	3
DNA dissolving reagent	7 (NR7025)	1
Suspending reagent	10 (NR10025)	1

Processing tube Sample tube: PT-8000 (8-hole, blue, 64pieces)

DNA tube: NT-8000 (8-hole, white, 64pieces)

KURABO INDUSTRIES LTD.

Bio-Medical Department

4-31, 2-chome, Kyutaro-machi, Chuo-ku

Osaka 541-8581, Japan

Phone: 81-6-6266-5010

Fax: 81-6-6266-5011

bio@ad.kurabo.co.jp