

DF-16

## Automated Plasmid DNA Extraction from *E. coli*

### Protocol

Collect the transformed *E. coli* into a 1.5 ml microtube, and pelletize

↓ ← RDP mix (RDP + EDP-01) \*1 : 100 µl

Vortex (No cell clumps should be visible after resuspension of the pellet)

Flash spin down

↓ ← ADP : 100 µl

Slowly mix by inverting the tube 5 times (Do not shake vigorously) \*2

Flash spin down (Do not leave the sample more than 5 min at this step)

↓ ← NDP : 140 µl

Slowly mix by inverting the tube 5 times (Do not shake vigorously) \*2

↻ ↓ 18,000 x g (14,100 rpm), 10 min, RT

Dispense 320 µl of LDP \*3 into a new 1.5 ml microtube

Transfer the supernatant (about 330 µl) to the 2 ml microtube\*4 with LDP

↓

Vortex (maximum speed) : 30 sec & Flash spin down

↓

Set into the device  
Protocol: PLASMID  
(Elution volume : 50 µl \*5)

\*Please refer to Quick Start Guide or operation manual  
to know how to set sample tube.

1. Apply the lysate into the cartridge
2. Pressurizing
3. Wash 2 times by Wash Buffer (WRT)
4. Add selected volume of Elution buffer and elute plasmid DNA into collection tube.

Plasmid DNA

\*1 Before starting an extraction experiment, add total amounts of EDP-01 to RDP bottle, and mix well. In the case of storing RDP mix, it is recommended to preserve it under refrigeration (2-8°C) and use within 6 months.

\*2 After addition of ADP or NDP, immediately mix by inverting the tube 5 times. Vigorous mixing results in the co-purification of much of genomic DNA. Too slow mixing causes inadequate blending of liquids, resulting in deterioration in the yield of plasmid DNA.

\*3 Add 44 ml of >99% ethanol into the bottle and mix well by gently inverting the bottle before use.

\*4 Following microtube are recommended.  
#BM4020  
(BM instrument co., ltd)  
#72.695.700,  
#72.695.500S  
(SARSTEDT)

\*5 The volume of the eluate from each cartridge is 100 µl. The volume of CRT can be reduced to 50 µl, but in that case, elution efficiency might be decreased.

## Results

---

The yield of plasmid DNA / Protein contamination : A260/280  
/ Chaotropic salt contamination : A260/230

Kit	A260/280	A260/230	Yield	S.D.
QuickGene	1.98	2.23	29.5 µg	4.07

N=4

Common protocol is usable for the following

---

Fosmid