

Specification

Product Name	Contents	Volume (for 48 tests)
GeneFields® Animal 5	PCR Oligo Mix	220μl × 1 Tube
	Coloring Buffer	1,600μl × 1 Tube
	DNA Strip	48 pieces × 1 Package
	User Manual & Interpretation Card	1 Set

Product Name	Contents	Volume (for 12 tests)
GeneFields® Animal 5 Plus	Box 1 (for PCR & Detection)	PCR Enzyme 80 μl × 1 Tube PCR Oligo Mix 60 μl × 1 Tube Coloring Buffer 410 μl × 1 Tube DNA Strip 12 pieces × 1 Package Interpretation Card 12 pieces × 1 Package
	Box 2 (for DNA extraction)	EDT 0.3 ml × 1 Tube MDT 3.1 ml × 1 Tube LDT 3.8 ml × 1 Tube WDT 20 ml × 1 Tube CDT 14 ml × 1 Tube Cartridges 12 pieces × 1 Package Waste Tubes 12 or 24 pieces × 1 Package Collection Tubes 12 pieces × 1 Package
	BioMasher II® ※1	12 pieces × 1 Package
	User Manual	1 booklet

- GeneFields® Animal 5 Plus Kit (**C** Type) (Code No.: GF-28-12**PC**) is supplied with Cartridges and Waste Tubes (24 pieces) to be used with centrifuge for DNA isolation.
- GeneFields® Animal 5 Plus Kit (**Q** Type) (Code No.: GF-28-12**PQ**) is supplied with Cartridges and Waste Tubes (12 pieces) to be used with KURABO QuickGene-Mini480 or QuickGene-Mini80 for DNA isolation.

※1. "BioMasher II®" is a registered trademark of Nippi. Inc., Tokyo, Japan.



KURABO INDUSTRIES LTD.

Bio-Medical Department, Advanced Technology Division

14-5, Shimokida-cho, Neyagawa-shi, Osaka 572-0823, Japan

Phone: +81-72-820-3079 FAX: +81-72-820-3095

www.kurabo.co.jp/bio/English/

June 2018