

Stock Preparation System for Laboratory

MONOPET

Easy Operation

The units already equip CPU and memory inside.
So just answer the question from MONOPET by pressing key.

Compact Body

The unit size small enough to set on the table.



Specification

Dosing Method	Liquid speed control, Gravimetric System
Balance	2,200/0.01g(±0.01g)
Dissolving Water	Hot water : Stainless Tank with temperature control (4 Litter) Cold water : Water Outlet
Speed	10cc/sec
Printer	Journal (24 char/line)
Dimensions	W405×D370×H740(mm)
Weight	50kg
Electric Power	AC100V±10V
Rated Power	1.5kW

Options

Stirring Unit (Option)	4 Bottles
------------------------	-----------



System Line up

AUCOLOR-TX10



Measuring the target colors, and making a recipe with selected dyestuffs or requested combination of dyestuffs.

AUKITCHEN-C100



Computer Color Kitchen for Laboratory. Dispense the dyestuffs based on the recipe. AUKITCHEN-C100 is good at dispensing small quantity.

AUKITCHEN mini



Weighing multiple dyestuffs simultaneously by using a dispenser. Mini volume measurement type with a turntable.

The appearance and specifications are subject to change for improvement without notice.

2019.11.



KURABO INDUSTRIES LTD. Advanced Technology Division

14-30, Shimokida-cho, Neyagawa, Osaka, 572-0823, Japan

TEL : +81-72-812-5205 FAX : +81-72-812-5261

URL : <https://www.kurabo.co.jp/el/world/en/> E-MAIL : dm_cis@kurabo-grp.com



Access here to read more!